Vesalius MCN *plus*

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**Quantum
Molecular
Resonance
Generator**

**Neurosurgery &
Spine Surgery**

Research and development
for your health

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NEUROSURGERY AND SPINAL SURGERY

Vesalius MCNplus is a modern mono and bipolar device, purposely developed for a wide type of operations. It finds its best application in procedures of Neurosurgery, SKull base and Spine Surgery.

Vesalius MCNplus is based on the Quantum Molecular Resonance Technology which creates high frequency fields combined in such a peculiar way in order to minimize the transfer of energy into a thermal form.

The quantized high frequency spectrum (4MHz with harmonics) generates quanta of energy with a sufficient value to break up the molecular bonds without increasing the kinetic energy, inducing only a modest temperature increase.

Characteristics

The CUT function, by using the QMR technology, induces breakage of molecular bonds while the temperature of the tissues barely increases. As a consequence, this function produces a precise cold cutting and causes a low lateral heat, therefore a minimal destruction of the tissue that surrounds the operation area, allowing an extremely delicate incision.

Activating the COAGULATION function, the MCNplus produces some waveforms off-resonance in a way to raise slightly the temperature.

As a consequence, instead of collapsing the vessels and necrotizing the area as the standard electro/radiosurgical units do, the effect of coagulation is obtained by triggering the denaturation process of the fibrinogen to form fibrin. The temperature of the target tissue increases more than when using the CUT function, but still remaining at modest level.

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Vesalius® is a patented product by Telea Electronic Engineering s.r.l.